



SEED STORAGE INCUBATOR



BSS Series

Shanghai Boxun Medical Biological Instrument Corp

Add: No. 599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

SEED STORAGE INCUBATOR >>

BSS Series

Product Features >

- The shell is made of cold-rolling steel with static spray plastics, the working chamber is made of mirror stainless steel, the shelves can be adjusted optionally;
- PMMA II control system, LCD screen displays temperature, time and humidity, over-temperature alarm function;
- Test hole is on both sides, internal diameter is 30mm;
- Centrifugal fan, brand compressor, fluorine free refrigerant;
- Equipped with mechanical lock,;
- Integrated type refrigeration system;
- Programmable design, it can be set 30 segments 99 cycles;
- Power supply recovery function, data will not be lost when power failure and system halted.



LCD Display



Shelf



Light



Model	BSS-150	BSS-250	BSS-400
Volume [L]	150	250	400
Power Supply	AC 220V+10V 50/60Hz		
Ambient Temperature [°C]	10~30		
Temp. Range without Humidifying [°C]	4~60		
Temp. Range with Humidifying [°C]	10~60		
Temp. Resolution [°C]	0.1		
Temp. Fluctuation [°C]	±0.5(heating); ±0.7(cooling)		
Temp. Uniformity at 25°C [°C]	±1		
Humidity Range (RH) [%]	20~60		
Power [W]	1300	1400	1750
Timer Range [min]	0~9999		
Internal Dimensions W×D×H [mm]	510*390*760	510*450*1100	600*640*1050
Housing Dimensions W×D×H [mm]	650*680*1400	650*740*1726	745*930*1695
Number of Shelves [pcs]	3	4	4
Load per Shelf [kg]	15	15	15
Net Weight /Gross Weight [kg]	114/131	132/162	162/194
Package Size [mm]	750*820*1590	750*880*1910	890*1120*1880