

## VERTICAL AUTOCLAVE



YXQ-A Series



YXQ-G Series

Shanghai Boxun Medical Biological Instrument Corp

Add: No. 599 Zhongqiang Road, Maogang Town,  
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

## VERTICAL AUTOCLAVE >>

### YXQ-A Series

#### Product Features >

- Made of SUS304;
- Hand wheel type quick-open cover;
- Microcomputer intelligent control;
- Self-expanding type seal ring;
- Exhaust cold air automatically;
- Inner circulation exhaust;
- Self-control when lack of water;
- Self-exhaust when over pressure;
- Pressure safety interlock;
- Over temperature protection;
- Buzzer reminds when sterilization ends;
- Equipped with G1/2A internal thread sample test hole

#### Three control modes:

- heating – sterilization - quick exhaust
- heating – sterilization - slow exhaust
- heating – sterilization - no exhaust



YXQ-A Series

### YXQ-G Series

#### Product Features >

- Drying function;
- Made of SUS304;
- Hand wheel type quick-open cover;
- Microcomputer intelligent control;
- Self-expanding type seal ring;
- Exhaust cold air automatically;
- Inner circulation exhaust;
- Self-control when lack of water;
- Self-exhaust when over pressure;
- Pressure safety interlock;
- Over temperature protection;
- Buzzer reminds when sterilization ends;
- Equipped with G1/2A internal thread sample test hole.

#### Work process of sterilization and drying:

- sterilization - natural cooling
- sterilization- automatic exhaust
- sterilization-exhaust-drain-drying



YXQ-G Series

Model	YXQ-50A	YXQ-70A	YXQ-100A
Volume [L]	50	75	100
Power [kW]	2.7	3.5	3.5
Drying Power [kW]	-	-	-
Power Supply [V, Hz]	220±10%, 50±2%		
Max. Working / Design Temp. [°C]	135 / 139		
Max. Working / Design Pressure [MPa]	0.22 / 0.25		
Timer Range [min]	4~120		
Internal Dimensions [mm]	Ø300×755	Ø400×720	Ø400×830
Housing Dimensions W×D×H [mm]	445×525×1120	485×605×1175	485×605×1290
Basket Dimensions [mm]	Ø280×210	Ø360×270	Ø360×230
Number of Baskets [pc]	3	2	3
Net weight /Gross weight [kg]	68/80	77/91	83/98
Package size [mm]	570×650×1240	610×730×1340	610×730×1490

Model	YXQ-50G	YXQ-75G	YXQ-100G
Volume [L]	50	75	100
Power [kW]	3.2	3.5	3.5
Drying Power [kW]	0.8	1	1
Power Supply [V, Hz]	220±10%, 50±2%		
Max. Working / Design Temp. [°C]	135 / 139		
Max. Working / Design Pressure [MPa]	0.22 / 0.25		
Timer Range [min]	4~120		
Internal Dimensions [mm]	Ø400×640	Ø400×720	Ø400×850
Housing Dimensions W×D×H [mm]	485×605×1125	485×605×1195	485×605×1290
Basket Dimensions [mm]	Ø360×230	Ø360×270	Ø360×230
Number of Baskets [pc]	2	2	3
Net weight /Gross weight [kg]	85/97	90/106	95/114
Package size [mm]	610×730×1290	610×730×1360	610×730×1490