



HORIZONTAL AUTOCLAVE



BXW-SD Series

Shanghai Boxun Medical Biological Instrument Corp

Add: No. 599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

HORIZONTAL AUTOCLAVE >>

BXW-SD Series

Product Features

- BRIGHT I intelligent operating system;
- 7-inch touch screen, real time display of parameters;
- With reservation function;
- Multi level permission management function;
- The system records device events and times on its own, and supports exporting and backing up USB data;
- Support over 1 million data records ;
- Built in 9 sterilization programs+3 melting and insulation programs;
- A 6-level cold air emission program can be established to ensure complete cold air emission Supports three-level exhaust setting ;
- Standard sample temperature probe ;
- Water level detection+anti dry burning low water level protection function;
- Equipped with steam generator and water pump, automatic water addition;
- Optional: printer, stainless steel casing.



Control Panel



Model	BXW-280SD-AI	BXW-360SD-AI	BXW-500SD-AI
Volume [L]	280	360	500
Door Structure	spoke	spoke	spoke
Power [kW]	12.5	15.5	17
Power Supply [V, Hz]	380, 50	380, 50	380, 50
Temp. Range [°C]	109~131	109~131	109~131
Internal Dimensions [mm]	Ø600×1050	Ø700×1050	Ø700×1450
Housing Dimensions W×D×H [mm]	1000×1515×1803	1000×1515×1803	1000×1915×1803
Net Weight [kg]	600	650	670
Gross Weight [kg]	650	700	720
Package Size [mm]	1600×1000×2000	1600×1000×2000	2000×1026×2000

Model	BXW-280SD-GI	BXW-360SD-GI	BXW-500SD-GI
Volume [L]	280	360	500
Door Structure	spoke	spoke	spoke
Power [kW]	13.5	16.5	18
Power Supply [V, Hz]	380, 50	380, 50	380, 50
Temp. Range [°C]	109~131	109~131	109~131
Internal Dimensions [mm]	Ø600×1050	Ø700×1050	Ø700×1450
Housing Dimensions W×D×H [mm]	1000×1515×1803	1000×1515×1803	1000×1915×1803
Net Weight [kg]	600	650	670
Gross Weight [kg]	650	700	720
Package Size [mm]	1600×1000×2000	1600×1000×2000	2000×1026×2000