

HIGH TEMPERATURE OVEN



BGG Series (400°C/500°C)

Shanghai Boxun Medical Biological Instrument Corp

Add: No. 599 Zhongqiang Road, Maogang Town,
Songjiang District Shanghai, China

Tel: 021-56980111 Fax: 021-56303876

E-mail: export@shbxyl.com

web: www.sh-bx.com



www.sh-bx.com

HIGH TEMPERATURE OVEN »

BGG Series

Product Feature

- The shell is made of cold-rolling steel with static spray plastics, the working chamber is made of mirror stainless steel, the shelves can be adjusted at will;
- Double press buckle hinge door lock makes the seal tighter;
- PMMA II control system, LCD display, timing, self-tuning function, accurate temperature control;
- Power supply recovery function, data will not be lost when power failure and system halted.



Control Panel



Hinged Door Lock

Model	BGG-38	BGG-78	BGG-148	BGG-248	BGG-39	BGG-79	BGG-149	BGG-249
Volume [L]	30	70	136	225	30	70	136	225
Temp. Range [°C]	RT+20~400	RT+20~400	RT+20~400	RT+20~400	RT+20~500	RT+20~500	RT+20~500	RT+20~500
Temp. Resolution [°C]	0.1	0.1	0.1	0.1	0.1	0.1	0.1	0.1
Temp. Fluctuation [°C]	±2	±2	±2	±2	±2	±2	±2	±2
Power [W]	1800	2050	2500	3000	2400	3400	4400	5500
Power Supply [V, Hz]	220, 50	220, 50	220, 50	220, 50	220, 50	220, 50	380, 50	380, 50
Timer Range [min]	0~9999	0~9999	0~9999	0~9999	0~9999	0~9999	0~9999	0~9999
Internal Dimensions W×D×H [mm]	300×300×350	400×400×450	450×550×550	500×600×750	300×300×350	400×400×450	450×550×550	500×600×750
Housing Dimensions W×D×H [mm]	670×480×700	770×580×800	820×730×900	870×780×1100	670×480×700	770×580×800	820×730×900	870×780×1100
Number of Shelves [pc]	2	2	2	2	2	2	2	2
Net Weight [kg]	61	86	108	132	61	86	108	132
Gross Weight [kg]	87	110	138	165	87	110	138	165
Package Size [mm]	800*660*910	900*763*1000	950*930*1100	1000*950*1300	800*660*910	900*763*1000	950*930*1100	1000*950*1300